

FILE 'CA' ENTERED AT 23:55:57 ON 25 APR 2008
USE IS SUBJECT TO THE TERMS OF YOUR STN CUSTOMER AGREEMENT.
PLEASE SEE "HELP USAGETERMS" FOR DETAILS.
COPYRIGHT (C) 2008 AMERICAN CHEMICAL SOCIETY (ACS)

FILE 'BIOSIS' ENTERED AT 23:55:57 ON 25 APR 2008
Copyright (c) 2008 The Thomson Corporation

FILE 'MEDLINE' ENTERED AT 23:55:57 ON 25 APR 2008

=> polyethylene glycol and (purified or purifying) and ion exchange
L1 412 POLYETHYLENE GLYCOL AND (PURIFIED OR PURIFYING) AND ION EXCHANGE

=> l1 and ion exchange chromatography
L2 235 L1 AND ION EXCHANGE CHROMATOGRAPHY

=> l2 and (branched polymer or multi-armed polymer)
L3 0 L2 AND (BRANCHED POLYMER OR MULTI-ARMED POLYMER)

=> l2 and purified polymer
L4 0 L2 AND PURIFIED POLYMER

=> (branched polymer or multi-armed polymer)
L5 3188 (BRANCHED POLYMER OR MULTI-ARMED POLYMER)

=> l5 and (purified or purifying)
L6 43 L5 AND (PURIFIED OR PURIFYING)

=> l6 and (polyethylene glycol or peg)
L7 1 L6 AND (POLYETHYLENE GLYCOL OR PEG)

=> d l7 bib abs

L7 ANSWER 1 OF 1 CA COPYRIGHT 2008 ACS on STN
AN 142:336828 CA
TI Method for preparing branched polyethylene glycol
IN Su, Zhiguo; He, Minglei
PA Institute of Process Engineering, Chinese Academy of Sciences, Peop. Rep.
China
SO Faming Zhuanli Shenqing Gongkai Shuomingshu, 13 pp.
CODEN: CNXXEV
DT Patent
LA Chinese
FAN.CNT 1

	PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
	-----	----	-----	-----	-----
PI	CN 1461762	A	20031217	CN 2002-120740	20020530
PRAI	CN 2002-120740		20020530		
AB	The method comprises reacting lysine (diamino acid, or polyamino acid) with HCl in anhydrous ethanol, esterifying to obtain lysine Et ester HCl; oxidizing methyl-polyethylene glycol with MnO2 at room temperature overnight then with 3% H2O2 for 24 h, separating on Bio-Rad Ag1*2 resin column with 0.02M HCl as eluent to obtain carboxymethylated mPEG; allowing to react lysine Et ester HCl (at a molar ratio of 2-4:1) in dichloromethane in the presence of triethylamine, dicyclohexylcarbodiimide, and N-hydroxysuccinimide at room temperature for 24 h, separating to obtain crude product; purifying on Biogel P100 5*50 column with water as eluent, extracting with dichloromethane, and recrystg. in ethanol.				

=>